







Specific Environment file








Specific Environment file





Other required software: Git



Required to get last development packages for the practicals

(Mac and Linux users may already have it on their computer!)
To check this enter « git --version » on your command line

(base) C:\Users\weber\Labo\Programmes Python\Git_others\FemtoUP2021\installation>git --version git version 2.26.0.windows.1





Creating and using the environment

In the anaconda prompt

Create a new environment

```
(base) C:\Users\weber\Labo\Programmes Python\Git_others\FemtoUP2021\installation>dir
Le volume dans le lecteur C n'a pas de nom.
Le numéro de série du volume est E425-BA9B
Répertoire de C:\Users\weber\Labo\Programmes Python\Git_others\FemtoUP2021\installation
                                     66 .gitattributes
                                 1 849 .gitignore
                            478 env_install.yml
268 env_install_pymodaq.yml
1 061 958 install PyMoDAQ Femto.odp
1 065 429 install.odp
             21:24
             22:16
                               759 809 install.pdf
                                   578 install_linux.pdf
                                   093 LICENSE
                                     73 README.md
                                       3 162 601 octets
                3 Rép(s) 301 847 945 216 octets libres
(base) C:\Users\weber\Labo\Programmes Python\Git_others\FemtoUP2021\installation>conda env create -f env_install_pymodaq.yml
```

(base) C:\Users\weber\Labo\Programmes Python\Git_others\FemtoUP2021\installation>conda env create -f env_install_pymodaq.yml

• Activate it —

conda activate pymodaq_femto

- Try the Retriever installation>retriever





Now is your turn!



